

Product Name : Algae Culture Agar	Product Code : NLAB-CHEMICALSIND21012
Description :	
Algae Culture Agar	
Technical Specification :	
Algae Culture Agar	
<p>Algae Culture Agar is recommended for the isolation and cultivation of algae from soil, water and sewage. Also for carrying stock cultures of algae used in the bioassay of algicidal chemicals. Algae Culture Agar is used for maintaining stock cultures of algae used in the bioassay of algaecide chemicals. The medium provides all necessary nutrients for good growth of algae but does not provide for other than minimal growth of bacteria. Stock cultures are prepared by inoculating the surface of slants with the algal culture and incubation at room temperature under a suitable light source.</p>	
Quality Control	
Appearance: Off-white to light yellow homogeneous free flowing powder	
pH: 6.80-7.20	
Gelling: Firm, comparable with 1.5% Agar gel	
Reaction: Reaction of 1.69% w/v aqueous solution at 25°C. pH : 7.0±0.2	
Colour and Clarity of prepared medium: White coloured clear to slightly opalescent gel forms in Petri plates.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003