

Product Name :
Alkali Aggregate Reaction Bath

Product Code :
NLAB-AGRETLAB160027



Description :

Alkali Aggregate Reaction Bath

Technical Specification :

The Alkali Aggregate Reaction Bath is used to keep 25 x 25 x 285 mm samples in NaOH (Sodium Hydroxide) or in any other solution at a specified temperature. The temperature can be adjusted between ambient to 80°C by using the digital controller with 2°C accuracy. The chamber is made of stainless steel and has a capacity for 36 specimens (25 x 25 x 285 mm).

Supplied

complete with 3pcs. stainless steel alkali specimens cans. Each can has a stainless steel hanger which can hold 3pcs. 25 x 25 x 285 mm specimens. Additional alkali specimens cans should be ordered separately.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003