

<b>Product Name :</b> Aluminium Calorimeter	<b>Product Code :</b> NLAB-EDUSCHOOLEQ1730189
<b>Description :</b>	
Aluminium Calorimeter	
<b>Technical Specification :</b>	
<p>This product is widely used for measuring the heat of chemical reactions or physical changes.</p> <p>Offered calorimeters are manufactured from optimum quality aluminum with 75 x 50 mm and 100 x 75 mm diameter.</p> <p>Furthermore, our product undergoes various quality assurance procedures, to ensure high standards of quality.</p> <p>Science Lab Supplies offer high quality Aluminum Calorimeter to our valuable customers.</p>	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003