

Product Name :
Amberlite CG-50 Hydrogen Form

Product Code :
NLAB-CHEMICALSIND11001

Description :

Amberlite CG-50 Hydrogen Form

Technical Specification :

Amberlite CG-50 Hydrogen Form

Grade: Laboratory Reagent (LR)

Capacity: 250gm, 1Kg

Description: Pale beige to white fine powder.

Exchange capacity: Min 10 meq/g

Particle size (100 – 200 mesh): Min 95%

Moisture About: 10.0%

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003