

Product Name :
Ambulance Stretcher Cum Wheelchair

Product Code :
NLABMEDICAL-M2110001



Description :

Ambulance Stretcher Cum Wheelchair

Technical Specification :

Ambulance Stretcher Cum Wheelchair

- * This Ambulance Stretcher can be converted into a chair.
- * The angle of the stretcher can be adjusted according to different requirements.
- * Suitable to transfer the patients in a limited place.
- * Such as lift in the hospital. ambulance, city road, etc.
- * Made of aluminum alloy material, characterized by its being anti-corrosion.
- * Use-safely and easy for the sterilization and cleaning.

Technical Parameter :

- * High Position (L x W x H) : 196 x 55 x 86cm

* Low Position : 196 x 55 x 25cm

* Packing Size (1 pc) : 198 x 62 x 26cm

* Max Angle of the Back : 85°

* Load Bearing : ?159kg

* N.W. (1 pc) : 34 kg

* G.W. : 40 kg

We are Exporters and Suppliers of Ambulance Stretcher cum Wheelchair at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003