

**Product Name :**  
Amplitude Modulation And Demodulation Kit

**Product Code :**  
NLAB-ENGINEERINGLB50004



**Description :**

Amplitude Modulation And Demodulation Kit

**Technical Specification :**

Amplitude Modulation And Demodulation Kit For Electrical Lab Training.

Amplitude modulation is implemented by means of a variable transconductance four quadrant multiplier designed around CA3086. On board carrier and modulating signals along with power supplies are provided. A detector and a low pass filter are provided to demodulate the signal. Modulation Index, Depth of Modulation, etc., are studied using this trainer.

**HARDWARE SPECIFICATION :**

- Carrier Generator 100 KHz
- Modulating signal generatorat 1KHz
- amplitude modulator
- amplitudedemodulator
- Inbuilt Power Supply
- Enclosed in an ABS plastic cabinet with cover
- set of 2 mmpatchchords manual

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003