

**Product Name :**  
Ampule Filling And Sealing Device( Manual)

**Product Code :**  
NLAB-SCIENTIFICEQ28003



**Description :**

Ampule Filling And Sealing Device( Manual)

**Technical Specification :**

The device has a Compact design with a cGMP model.

The Ampoule gas is centring for the period of gas flushing and filling.

There is Possibility of gas flushing prior to, during and post filling.

AISI 316L dosing syringes.

“No ampoule”. It does not need any filling function.

The machine is equipped with pre-warming station.

Input Specification:15 to 30 Ampoules per minute

Container Dia.: 10mm to 21.5mm Container Height: 60mm to 140mm

Filling Range :1 ml to 20 ml

Filling Accuracy :A± 1% from filling volume Aluminum Profile Cabinet covering full machine

Power Specification: 1 H.P., 3 Phase, 4 Wire System, 50 Hertz.

Nitrogen Pressure: 1kg/Cm2 Oxygen Pressure Required: 1kg/Cm2 LPG Pressure Required: 0.5 kg/Cm2

Gas for Sealing: Oxygen + LPG

Optional Accessories: Programmable Logic Control (System) with MMI Aluminum Profile Cabinet covering full machine

Overall Dimension:915mm (L) X 610mm (W) X 965mm (H)

Electrical characteristics can be adjusted as per customer requirements.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003