

Product Name :
Anaesthesia Machine Diamond Economy

Product Code :
NLABMEDICAL-M2230001



Description :

Anaesthesia Machine Diamond Economy

Technical Specification :

Anaesthesia Machine Diamond Economy

The main body is made of High-strength engineering plastic rack, light, beautiful and corrosion resistant

It can be applied to both Adult & Child

The gas source is O₂: 0.28~0.6 MPa

Flowmeter of O₂: 0.05~1.0L/min is 1.1~10L/min

The flow rate of rapid oxygen supply is 25~75L/min

Vaporizer Among Halothane, Enflurane, Isoflurane, and Sevoflurane, anyone can be chosen for application as required by the customer.

The respiratory circuit has working mode: all-close, semi-closed, semi-open APL:~12.5kPa

Respiratory bellows for adults, bellows for children are supplied with the Machine

The tidal volume range is 0~1500 ml

We are Exporters and Suppliers of Anaesthesia Machine Diamond Economy at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003