

Product Name :
Anaesthesia Machine Major S

Product Code :
NLABMEDICAL-M2130001



Description :

Anaesthesia Machine Major S

Technical Specification :

Capacious drawer for keeping anaesthesia kit/accessories

Extended platform at rear for two 10ltrs. (Water Cap.) Cylinders

Large roller bearings, antistatic castors, for easy mobility and brakes

Top tray provided at eye level for monitoring equipments

Hook for breathing circuit

Clean stainless steel top

Provision to fit ventilator

Emergency drug tray fitted at convenient place
Unitised rota meter block for oxygen & nitrous oxide with fine adjustment controls & colour coated knob and luminous back plate. Oxygen flow meter calibrated from 100cc/min to 10 ltrs./ Minute & Nitrous oxide flow meter from 200cc/mm to 12 liters/min

One gold mantype vaporizer for halothene, provision to mount three vaporizer at a time.

Non return cum blow off valve in circuit ensures that pressure does not build-up beyond 200cm. Water column & minimises risk of back flow of gases.

60 psl outlet for driving ventilator

Change over system open circuit or close circuit

Emergency Oxygen flush valve at table level to give 35-60 ltrs/min for emergencies

Cylinder yokes: pin indexed yokes two each of oxygen & nitrous oxide

63 mm diameter, colour coded gauges fitted at convenient angle for better visibility & clarity, one gauge for O2 & one gauge for

N2O shown the line pressure

High efficiency regulators for reliability, two each for O2 & N2O

One Maggills or Brain Circuit

One B.P Apparatus.

Safety Features

Hypoxic guard

Emergency oxygen flush

Alarm to provide audible indication of O2 supply Pressure failure

Nitrous oxide shut off in the event of O2 failure

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003