

**Product Name :**  
Analog Tensile Testing Machine

**Product Code :**  
NLAB-ENGINEERINGLB31007



**Description :**

Analog Tensile Testing Machine

**Technical Specification :**

This machine is used to held various tests like Compression, Cross breaking / Bending, Shear, Linear Bonding / Gripping / Peeling strength of adhesives & tapes etc.

It can test materials like rubber, wood, plastic, rexins, teflon, cardboard, composites, ferrous metals, non-ferrous metals etc.

Wide range of Grips, Attachments & Special fixtures are available to suit varieties of material and shapes. Continuous roll autographic Load Vs Elongation recorder.

Continuous roll autographic Load Vs Elongation recorder. Elongation scale with least count 1 mm.

Elongation scale with least count 1 mm. Three measuring ranges for accurate testing.

Three measuring ranges for accurate testing. A number of straining speeds are available for proper selection.

A number of straining speeds are available for proper selection. Machine with AC or DC drive motor is available.

Machine with AC or DC drive motor is available. Various Models – Maximum Capacity from 200 N to 50 kN.

Various Models – Maximum Capacity from 200 N to 50 kN. Loading accuracy well within  $\pm 1\%$

Loading accuracy well within  $\pm 1\%$  Hand operated cutting press with different sizes dies are available for cutting dumbbell specimen.

Hand operated cutting press with different sizes dies are available for cutting dumbbell specimen.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003