


Product Name :
Analog to Digital Converter Kit

Product Code :
NLAB-SCHOOLABS13010



Description :

Analog to Digital Converter Kit

Technical Specification :

Shows simple principles of analog signal conversion. Continuously converted to a 3-bit digital signal, each bit is displayed by a 3-color LED. The results are clearly readable and easy to understand.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003