

Product Name :
Analogus to Digital Converter

Product Code :
NLAB-VOCATIONALTVET14003



Description :

Analogus to Digital Converter

Technical Specification :

Analogus to Digital Converter.

L7805 + 5V D.C. at 250mA, IC Regulated power supply internally connecated.

- 5V D.C. at 50mA,

IC regulated power supply internally connected .

+ 10V D.C. at 50mA ,IC Regulated power supply internally connected.

0-15V D.C. at 50mA, continusly variable power supply .

D.C. Voltmeter , 65mm continusly dial to read 0-15V.

Voltage comparator Quad 2 -input NOR gate.

Qual 2-input AND gate. Qual 4-input EX-OR gate.

BCD-to-7 Segment Decoder/Driver.

Dual J-K Flip-Flop.

Bit Binary Ripple Counter. Timer.

Continuous monitoring of analog signals on a voltmeter and of digital on 7 segment display.

Seven LEDs to display the of various important points.

Electrical Work Bench Lab Analogus to Digital Converter Manufacturer, Supplier and Analogus to Digital Converter Exporter in India

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003