

**Product Name :**  
Analogus to Digital Converter

**Product Code :**  
NLAB-VOCATIONALTVET14003



**Description :**

Analogus to Digital Converter

**Technical Specification :**

Analogus to Digital Converter.

LTB-812 + 5V D.C. at 250mA, IC Regulated power supply internally connecated.

- 5V D.C. at 50mA,

IC regulated power supply internally connected .

+ 10V D.C. at 50mA ,IC Regulated power supply internally connected.

0-15V D.C. at 50mA, continusly variable power supply .

D.C. Volmeter , 65mm continusly dial to read 0-15V.

Voltage comparator Quad 2 -input NOR gate.

Qual 2-input AND gate. Qual 4-input EX-OR gate.

BCD-to-7 Segment Decoder/Driver.

Dual J-K Flip-Flop.

Bit Binary Ripple Counter. Timer.

Continuous monitoring of analog signals on a voltmeter and of digital on 7 segment display.

Seven LEDs to display the of various important points.

Electrical Work Bench Lab Analogus to Digital Converter Manufacturer, Supplier and Analogus to Digital Converter Exporter in India

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003