


**Product Name :**  
Analytical Balance

**Product Code :**  
NLAB-SCIENTIFICEQ14002



**Description :**

Analytical Balance

**Technical Specification :**

**Large glass drought shield** with 3 sliding doors for easy access to the items being weighed

- Large LCD display, • stainless steel, Ø 80 mm pan

**Automatic internal adjustment**

**Data interface RS-232:** To connect the balance to a printer, PC or network

**Dosage aid:** High-stability mode and other filter settings can be selected

Tare/print function. mode of selection in gm/kg/pound. 4 digits, printed on calibration protocol freely programmable Permissible ambient temperature 10 °C/30 °C Sensitivity available:

- 0.1gm
- 0.01gm

- 0.001gm
- 0.0001gm/0.1mg

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003