


Product Name :
Analytical Balance

Product Code :
NLAB-SCIENTIFICEQ14002



Description :

Analytical Balance

Technical Specification :

Large glass drought shield with 3 sliding doors for easy access to the items being weighed

- Large LCD display, • stainless steel, Ø 80 mm pan

Automatic internal adjustment

Data interface RS-232: To connect the balance to a printer, PC or network

Dosage aid: High-stability mode and other filter settings can be selected

Tare/print function. mode of selection in gm/kg/pound. 4 digits, printed on calibration protocol freely programmable Permissible ambient temperature 10 °C/30 °C Sensitivity available:

- 0.1gm
- 0.01gm

- 0.001gm
- 0.0001gm/0.1mg

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003