

**Product Name :**  
Analytical Balance

**Product Code :**  
NLAB-SCIENTIFICEQ14003



**Description :**

Analytical Balance

**Technical Specification :**

- A precision balance suitable for quick and accurate weighing. It has a short armed triangular type beam with zero at the center and 50 graduated notches on either side of zero for rider use.
- It is highly suitable for research work in universities and industrial laboratories specially designed to maintain constant sensitivity under varying load and acute self-acting parallel arrestment of the beam.
- Complete balance housed in a wooden case with leveling screws and clear glass front and access from both the sides.
- Capacity : 200g
- Sensitivity : 0.1 mg
- Pan: Stainless steel pans of 75mm diameter, concave upwards.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003