

Product Name :
Analytical Balance

Product Code :
NLAB-SCIENTIFICEQ14003



Description :

Analytical Balance

Technical Specification :

- A precision balance suitable for quick and accurate weighing. It has a short armed triangular type beam with zero at the center and 50 graduated notches on either side of zero for rider use.
- It is highly suitable for research work in universities and industrial laboratories specially designed to maintain constant sensitivity under varying load and acute self-acting parallel arrestment of the beam.
- Complete balance housed in a wooden case with leveling screws and clear glass front and access from both the sides.
- Capacity : 200g
- Sensitivity : 0.1 mg
- Pan: Stainless steel pans of 75mm diameter, concave upwards.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003