

Product Name :
Analytical Balances with External Calibration

Product Code :
NLABMEDICAL-M2960001



Description :

Analytical Balances with External Calibration

Technical Specification :

Analytical Balance is also known as Laboratory Digital Balance

Automatic external calibration is available in this model

It is LCD Display with backlight

The Standard RS 232 C Interface is available in this model

This is based on E.M.F.C. Technology

The balance has S.S weighing pan, Level Indicator, adjustable fit

Full tare range up to a maximum capacity of the Balance

Repeatability : ± 0.1 mg

Linearity : ± 0.2 mg

Pan size: 85 \varnothing mm

Stabilization Time: 3-4 Sec

Communication: RS 232 Interface

Calibration : External calibration

Weight Units : Multi unit

Power : 230 V / 11V AC or 120V/ 11V AC

Display : LCD Display with back light

Working Temperature : 10° C – 40° C

We are Manufacturer, Exporters and Suppliers of Analytical Balances with External Calibration at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003