

Product Name :
Anchors Horizontal Lifeline System

Product Code :
NLAB-SAFETYLABP20021



Description :

Anchors Horizontal Lifeline System

Technical Specification :

Webbing : Fray-proof, dope-dyed polyester webbing

Strength : 15 kN

Lifeline Length : 20 metres

Ratchet : Steel, zinc coated, rated to 50 kN

Features

Equipped with a Ratchet Tensioner for quick and easy tensioning of the lifeline between two supports

Fall arrest rated for two users, provides job site flexibility

Built in shock absorber reduces anchorage loadings to preserve the structure

Complete kit includes retaining bag and two tie-off anchor slings

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003