

**Product Name :**  
Anchors Monopod with Winch

**Product Code :**  
NLAB-SAFETYLABP20013



**Description :**

Anchors Monopod with Winch

**Technical Specification :**

Specially developed for stationary use in safety and rescue cases for two persons as well as for lifting and lowering loads.

Removable extension arm, adjustable height.

Suitable for fast relocation.

Stable thanks to solid construction, high quality, Height 1.8 m, cantilever reach 1.25 m.

Safe Working Load Limit up to 300 kg with max 20m rope length.

V2A stainless steel design.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

