

Product Name :
Anchors Polyester Sling

Product Code :
NLAB-SAFETYLABP20018



Description :

Anchors Polyester Sling

Technical Specification :

Tie of Anchor slings made of polyester

Strength: 22 kN.

Length: 1m/1.5m/2 m.

Round Sling with D-ring at one end.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003