

Product Name :
Anchors SS Wire Rope Sling

Product Code :
NLAB-SAFETYLABP20020



Description :

Anchors SS Wire Rope Sling

Technical Specification :

Made of 8 mm Stainless Steel grade 304

Available Length: 1 m/1.5 m/1.8 m/2 m or custom length.

Loops at both ends.

Strength: 15 kN.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003