

Product Name :
Anchors Telescopic Anchor Device

Product Code :
NLAB-SAFETYLABP20011



Description :

Anchors Telescopic Anchor Device

Technical Specification :

Opening Gate: 90 mm

Breaking Strength: 25 kN

Weight: Light weight

Anchor Device: Aluminium Alloy/Steel

Length: Minimum 2 Mtr.; Maximum 6 Mtr.

Maximum Fixing Height: 7.6 Mtr. (including human height)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

