

Product Name : Anchors Tripod	Product Code : NLAB-SAFETYLABP20014
---	---



Description :
Anchors Tripod
Technical Specification :
Auto-locking system at the head to ensure stability of the legs.
Built-in anchor point : Eye bolt offers secondary anchor point for fall arrest system
Material: High-strength Aluminum Alloy.
Adjustable telescoping legs : Easy pin style adjustment allowing precise height setting.

Portable Device for Confined Space - Entry and Retrieval

Strength: 22 kN.

Safe Working Load Limit 200 kgf.

Anti-skid rubber foot-pads and safety straps prevent legs from spreading apart.

Height adjustments: From 1.45m to 2.15m.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003