

**Product Name :**  
Anchors Webbing Sling

**Product Code :**  
NLAB-SAFETYLABP20019



**Description :**

Anchors Webbing Sling

**Technical Specification :**

Anchor slings made of 44 mm polyester webbing

Strength: 22 kN.

Length: 1m/1.5m/2 m.

D-rings at both end.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003