


Product Name :
Aneroid Barometer

Product Code :
Weather0002



Description :

Aneroid Barometer

Technical Specification :

Plastic base approx. 160mm Dia.

100mm dial graduated 28 to 31 ϵ^3 of mercury and 960 to 1060 millibars with plated bezel, cover glass and index

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-0171-2643080, 2601773