

Product Name :
Anti Static Full Body Harnesses

Product Code :
NLAB-SAFETYLABP20025



Description :

Anti Static Full Body Harnesses

Technical Specification :

Materials: Anti-static webbing to prevent the risk of build up charge.

Stitching: High-strength polyester thread stitching

Anchor point: One dorsal D-ring and Two textile chest attachment loop.

Special Features: Ideal use for explosive atmosphere when the risk of an electrostatic build up and sudden discharge, igniting.

Adjustability: Size adjustment strap only at chest & thighs.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003