

Product Name : Anti-Static Rebreathing Bag	Product Code : NLABMEDICAL-M3250018
--	---

--

Description :
Anti-Static Rebreathing Bag India, Anti-Static Rebreathing Bag manufacturer India, suppliers, Medical Anti-Static Rebreathing Bag India, Hospital Anti-Static Rebreathing Bag India

Technical Specification :
Anti-Static Re-Breathing Bag, Black Color, Size: ½ ltr.
Anti-Static Re-Breathing Bag, Black Color, Size: 1 ltr.
Anti-Static Re-Breathing Bag, Black Color, Size: 1.5 ltr.
Anti-Static Re-Breathing Bag, Black Color, Size: 2 ltr.
Anti-Static Re-Breathing Bag, Black Color, Size: 3 ltr.

Naugralabequipments
Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003