

Product Name :
Apparatus For Methylene Blue Test

Product Code :
NLAB-EDUCATIONTVET13021



Description :

Apparatus For Methylene Blue Test

Technical Specification :

The set includes:

Complete agitator with blades, diameter 70 mm, 400/700 rpm, 220 V / 50 Hz

Tripod for agitator

Filter package, diameter 90 mm

Glass rod, 300 mm

Container, 2000 ml

Plastic tray 200x150x80 mm

Methylene blue, 100 g

Kaolinite, 500 g

Burette 50 ml x 0.1 ml

Burette holder with clamp

Burette funnel

Amber glass reagent bottle, 500 ml

Liquid dispenser, 0.5 to 5.0 ml

Plastic cap

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003