

<b>Product Name :</b> Ascarite AR 1.6-3 mm Acid Washed Absorption	<b>Product Code :</b> NLAB-CHEMICALSIND13001
<b>Description :</b> Ascarite AR 1.6-3 mm Acid Washed Absorption	
<b>Technical Specification :</b> Capacity:100 gm, 250 gm Grade:Laboratory Reagent (LR) Capacity for CO2	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003