

Product Name :
Asphalt Coring Machine

Product Code :
NLAB-ENGINEERINGLB31008



Description :

Asphalt Coring Machine

Technical Specification :

The Laboratory Asphalt Coring Machine is capable of achieving a large range of spindle speeds to provide optimum performance when preparing test specimens, regardless of the properties of the material being sampled. All accessories are included (water swivel, drip pan, sample holder, and core barrels for coring 100 mm (4 inch) specimens. Other standard or custom barrels are also available.

The main feature that sets this coring machine above any other is the force/stroke control to optimize feeding rates for superior performance when coring material, which is not homogeneous.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003