

**Product Name :**  
Aspirator Bottles

**Product Code :**  
NLAB-SCIENTIFICEQ37002



**Description :**

Aspirator Bottles

**Technical Specification :**

Are made of Polypropylene and are therefore much lighter than glass. These are autoclavable, durable & dependable. The leakproof spigot screwed to the bottle body requires only a half turn to open or close. These bottles are used for storing or dispensing distilled water & other laboratory solutions. Threads should be totally disengage before autoclaving.

**Material :** Polypropylene

**Catalogue No .**

NTLR7.01

NTPt7.02

NTPt7.03

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003