

**Product Name :**  
Auditory Trainers

**Product Code :**  
NLABMEDICAL-M3220004

**Description :**

Auditory Trainers India, Auditory Trainers manufacturer India, suppliers, Medical Auditory Trainers India, Hospital Auditory Trainers India

**Technical Specification :**

Amplifier: IC-based output stage coupled with 2 transistor low-noise pre-amp-stage.

Freq. Range: 150 Hertz to 8 Khz.

Gain Control: Continuously variable dB Attenuator Max. Volume 130 dB SPL.

Tone Control: 2 Tone per channel: Hi-Low (Fixed)

Indicators: Separate indicators(VU Meters)

Switches: For 'Mike' and phono(induction loop system). For battery check.

Microphones: Electret-built in and external

Headphones: Dynamic, low impedance, Adjustable head-band and circumaural ear cushion.

Supply Voltage: 6V; 1.5x4 R20(large size) cells if on battery. 230V; 50Hz. A.C. mains with built in power supply.

Current Drain: 20 mA, Quiescent.

Max. Power: 700 mW, per channel

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003