

**Product Name :**  
Autoclave Bags

**Product Code :**  
NLABMEDICAL-M2280001



**Description :**

Autoclave Bags

**Technical Specification :**

Autoclave Bags

Autoclave Bags are made of medical-grade HDPE / PP

These are used contaminated medical waste disposal in sterilizers and incinerators. These bags have the highest melting point.

Autoclave Bags will not melt, even in superheated vacuum sterilizers reaching up to 285°F/141°C.

Designed for sterilization of sensitive products and components

Autoclave Bags can be printed with international biohazardous symbols and warning languages

Autoclave Bags meet the strength requirements of 165-gram dart drop standard as performed

These Bags are perfect for steam, Eto and gamma sterilization processes

The Colour available are Red, Yellow, Clear, plain clear with and without printing

We are Manufacturers, Exporters, and Suppliers of Autoclave Bags at tender and export competitive prices from India.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003