

Product Name :
Autoclave Pressure Cooker Aluminium Electric

Product Code :
NLABMEDICAL-M2240001



Description :

Autoclave Pressure Cooker Aluminium Electric

Technical Specification :

Autoclave Pressure Cooker is also known as Pressure Steam Sterilizers, Steam Autoclave, Portable Autoclaves Sterilizers

Autoclave Pressure Cooker is a part of Steam Sterilizing Equipment

Autoclave Pressure Cooker is Electrical and has a mirror finish

The body and lid have Seamless Construction from Extra Thick Aluminium Sheet Body & Lid for more safety and life long performance.

The Anti-Bulging & Fuel Saver Base is provided for Additional Strength, Fuel Saving and Time Saving.

The depression at the centre of the lid helps in Additional Strength.

The Gasket is ISI Marked food grade synthetic rubber gasket for better reliability, hygiene & longer life.

The Gasket Release System (GRS) & fusible plug type safety valve is available for Triple-Safety.

The pressure gauge is Colour-coded Pressure Gauge to indicate the safe and danger zones.

The repeating type of safety valve is provided for better Safety and longer & error-free performance.

The vacuum release & steam release valves are made of Stainless Steel for longer life and to avoid toxic deposits on valves due to oxidation.

The heating element is ISI Marked which is made of thick and heavy copper pipe with very thick chrome plating.

The Electric Cable has heavy electric wires to comfortably bear the heating load.

The 3 Pin Plug & Element Connector is of High quality for electrical safety.

The handles are made from heat resistant, impact tolerant, finest quality bakelite, fitted with Stainless Steel screws to avoid easy loosening of handles.

The inside tripod stand is made of Stainless Steel with extra height to avoid wetting of instruments.

The outer stand is well designed to provide proper stability to the Autoclave Pressure Cooker.

We are Manufacturer, Exporters and Suppliers of Autoclave Pressure Cooker Aluminium Electric at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003