

Product Name :
Autoclave Vertical For Laboratory

Product Code :
NLAB-LABORATORYEM29001



Description :

Autoclave Vertical For Laboratory

Technical Specification :

For sterilization of culture media, glassware, utensils instruments etc. in steam under pressure.

Steam lined construction.

Outer body is made of thick mild steel sheet finished with cream enamel painted.

Inner chamber is made of thick stainless steel.

Lid is made of mild steel lined with stainless steel & provided with a radial locking systems, worked by a paddle at the .

Neoprene rubber, hydraulically tested up to 40 psi.

The chamber is absolutely leak proof & can be operated at any selected point in between 5 to 20 pound per sq. inch.

It is provide with safety valve pressure gauge, pressure release valve, water level indicator with cord & plug.

Fitted with automatic cut off device for element for general safety & also fitted with indfoss piezostat.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003