

**Product Name :**  
Autoclave Vertical Semi Automatic

**Product Code :**  
NLAB-LABORATORYEM29006



**Description :**

Autoclave Vertical Semi Automatic

**Technical Specification :**

For Dry Sterilization of Autoclave Vertical Semi Automatic. Digital display of temperature and timer. Inner chamber, outer wall, boiler & steam jacket are made of S. S. 304 grade. Middle Chamber : Stainless Steel 304 Lid : Stainless steel 304 Radial Locking system : Made of steel bars Accessories included:

Automatic pressure control switch

Automatic Low water cut off arrangement

Vacuum Breaker

Pressure Gauge

Water level indicator

Safety Valve

Steam Release Valve

Paddle lifting Device

S. S. Basket

Special Features:-

Single panel feather touch LED microprocessor controller

Simultaneous display of temperature & time.

Timer starts automatically as soon as the set temperature is attained & stop the cycle automatically as soon the time is over.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003