

**Product Name :**  
Autoclave

**Product Code :**  
NLAB-SCIENTIFICEQ12003



**Description :**

Autoclave

**Technical Specification :**

Capacity: 200 Liter

Chamber Size(mm): 750\*1000

Chamber Material: SUS304

Designed Pressure: 0.24 Mpa

Designed Temp: 138°C

Working Pressure: 0.217 MPa

Working temp.: 50~134°C

Temp. Display accuracy: ±0.1°C

Timing range: 0~99h

Dynamic Pulse: Exhaust Times: 0~9 times; Exhaust

Temp.: 105~134°C

Dissolution Temp.: Range 60~100? (programmable); Temp. Holding Range 40~60°C (programmable)

SS Basket: 4pcs

Gross Weight: 250 kg

Power Connection: 380V, 50/60 Hz

Provide the warranty period of 24 months

Wet sterilization of micro-lab testing apparatus, Media and Buffer Solutions.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003