

Product Name :
Autoclave

Product Code :
NLAB-SCIENTIFICEQ12012



Description :

Autoclave

Technical Specification :

Average cycle time: 30 min

Chamber size (Volume/nominal): 80/75 Liter

Max psi: 37 psi

Chamber construction: Stainless Steel

Display: 4 digit LED

Pressure: 121°C at 1.5 bar

Power Connection: 240 V, 50/60 Hz

Provide the warranty period of 24 months

Sterilization of media and tools before use.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003