

Product Name :
Autoloader Collapsible Stretcher (Ambulance Stretcher)

Product Code :
NLABMEDICAL-M2110002



Description :

Autoloader Collapsible Stretcher (Ambulance Stretcher)

Technical Specification :

- * This automatic loading stretcher is our latest product, it's made of high-strength aluminum alloy thick-wall pipes.
- * The bed surface is made of foamed cushion and the adjustable back making the wounded feel comfortable.
- * The mechanical foldaway legs are controlled by the handles on both ends of the stretcher.
- * And the legs can be folded under the stretcher automatically when loading into the vehicle.
- * Also the stretcher can be locked by the fixing devices when on the ambulance car.
- * This type of stretcher is characterized by its being light.
- * Durable antiseptic easy for sterilization.

Technical Parameter :

- * High Position (L x W x H) : 190 x 55 x 92cm.

- * Low Position : 190 x 55 x 25cm.
- * Packing Size (1 pc) : 195 x 62 x 25.5cm
- * Max Angle of the Back : 75°
- * Load Bearing : ? 159kg
- * N.W. (1 pc) : 40 kg
- * G.W. : 46 kg

We are Exporters and Suppliers of Autoloader Collapsible Stretcher (Ambulance Stretcher) at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003