

**Product Name :**  
Autoloader Collapsible Stretcher (Ambulance Stretcher)

**Product Code :**  
NLABMEDICAL-M2110002



**Description :**

Autoloader Collapsible Stretcher (Ambulance Stretcher)

**Technical Specification :**

- \* This automatic loading stretcher is our latest product, it's made of high-strength aluminum alloy thick-wall pipes.
- \* The bed surface is made of foamed cushion and the adjustable back making the wounded feel comfortable.
- \* The mechanical foldaway legs are controlled by the handles on both ends of the stretcher.
- \* And the legs can be folded under the stretcher automatically when loading into the vehicle.
- \* Also the stretcher can be locked by the fixing devices when on the ambulance car.
- \* This type of stretcher is characterized by its being light.
- \* Durable antiseptic easy for sterilization.

**Technical Parameter :**

- \* High Position (L x W x H) : 190 x 55 x 92cm.

\* Low Position : 190 x 55 x 25cm.

\* Packing Size (1 pc) : 195 x 62 x 25.5cm

\* Max Angle of the Back : 75°

\* Load Bearing : ? 159kg

\* N.W. (1 pc) : 40 kg

\* G.W. : 46 kg

We are Exporters and Suppliers of Autoloader Collapsible Stretcher (Ambulance Stretcher) at tender and export competitive prices from India.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003