

Product Name :
Automatic Compression Testing Machine

Product Code :
NLAB-TESTINGLABM31003



Description :

Automatic Compression Testing Machine

Technical Specification :

The Automatic Compression Testing Machine has all features of Micro Controller based Compression Testing Machine. In addition, it has

The machine is available in 500kN, 1000kN, 2000kN, 3000kN ranges.

An automatic load release facility.

A facility for transmission of data to the serial port of a computer placed within 200 meters of the loading unit.

Support
for data transmission over distances up to several kilometers. This is useful for testing in the field. (Optional feature)

Measurements of Strain with a resolution of up to 0.001 mm (optional feature).

Measurement of Flexural Strength (optional feature)

Facility to enter the required pace rate, at the beginning of the run.

The load resolution is for is 0.1kN for machines up to 1000kN capacity & 1.0kN for the 2000kN & 3000kN

machines.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003