

Product Name : Automatic Dental Amalgamator	Product Code : NLABMEDICAL-M32100024
Description :	
Automatic Dental Amalgamator India, Automatic Dental Amalgamator manufacturer India, suppliers, Medical Automatic Dental Amalgamator India, Hospital Automatic Dental Amalgamator India	
Technical Specification :	
Automatic amalgamation with timer is a compact and functionally designed mixer which assures efficient and reliable triturating with a simple operation	
Excellent durability and mixing power with 4 capsules	
Dimensions: 13 cm. Width, 22cm. Length, 12 cm. Height, 3.5 kg Weight	
Power 220/ 230 Volts AC Consumption: 45 W	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003