

<b>Product Name :</b> Automatic Dental Amalgamator	<b>Product Code :</b> NLABMEDICAL-M32100024
<b>Description :</b> Automatic Dental Amalgamator India, Automatic Dental Amalgamator manufacturer India, suppliers, Medical Automatic Dental Amalgamator India, Hospital Automatic Dental Amalgamator India	
<b>Technical Specification :</b> Automatic amalgamation with timer is a compact and functionally designed mixer which assures efficient and reliable triturating with a simple operation  Excellent durability and mixing power with 4 capsules  Dimensions: 13 cm. Width, 22cm. Length, 12 cm. Height, 3.5 kg Weight  Power 220/ 230 Volts AC Consumption: 45 W	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003