

Product Name :
Automatic Digital Compression Testing Machine

Product Code :
NLAB -CONCRLAB120009



Description :

Automatic Digital Compression Testing Machine

Technical Specification :

Automatic test cycle with standard rapid approach

Graphical – LCD data acquisition and control system

Menu driven software for easy operation

Load vs time plot and instantaneous load rate displayed

Accurate load rate control within ± 5 from 1 kN/sec to 20 kN /sec

Class 1 starting from % 10 of the full range

Ready stress calculation for the standards of EN 12390-3, EN 12390-5, EN 12390-6, En 1338, EN 1340, EN 196

SI, Metric, Imperial units

Data storage up to 256 tests

RS 232 output for printer or PC

Free of charge PC software for test control and advanced report printout

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003