

**Product Name :**  
Automatic Flexure Tension Machine

**Product Code :**  
NLAB-CEMETLAB130038



**Description :**

Automatic Flexure Tension Machine

**Technical Specification :**

Wattage: 40 W

Dimensions: 510x1050x720 mm

Weight approx.: net 52 kg

Constant rate of loading

Dual range

Automatic stop at failure

Accepts flexural or briquette specimens

Single lever loading assures high accuracy

Automatic flexure/tension machine for cement specimens. 5000 N cap. 230 V, 50 Hz, 1 ph.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003