

Product Name :
Automatic Knife Sharpner (Shadon Type)

Product Code :
NLAB-LABORATORYEM28004



Description :

Automatic Knife Sharpner (Shadon Type)

Technical Specification :

Automatic Knife Sharpener (Shandon Type)

Key Features :-

SCM Program control, easy operation.

Multi purpose knife holder.

Automatic disc leveling.

Buffer adjusting mechanism for turning knife.

Convenient loading and unloading for knife slicing.

Low noise operation.

Technical Specifications :-

Knife Length : 110-250mm

Thickness : Less than 13mm

Grinding Angle : 22-45 deg C (half angle 11-22.5 deg C)

Timing Range : 1-60 min

Knife interval : 5,10,25 s and unilateral

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003