

Product Name :
Automatic Razor Sharpener (Arthur Thomas Type)

Product Code :
NLAB-ENGINEERINGLB95002



Description :

Automatic Razor Sharpener (Arthur Thomas Type)

Technical Specification :

Razor Sharpener is Arthur Thomas type heavy duty sharpening machine for razors upto 300 mm length.

It is equipped with large round hard surface plate with one side ground and other side polished for various sharpening operations.

Completed with universal knife holder having automatic reversing device and dust cover without abrasives.

Optional:-

Micro abrasives for coarse & fine honing are available in 200 gm

Packing at extra cost

Suitable for all models of sharpener

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments