

Product Name :
Automatic Razor Sharpener (Spencer Type)

Product Code :
NLAB-ENGINEERINGLB95001



Description :

Automatic Razor Sharpener (Spencer Type)

Technical Specification :

Automatic Razor Sharpener is a spencer type streamlined automatic microtome knife sharpener with electric timer coupled with simple operation for razors up-to 185 mm on both the wedges. Working mechanism encased in aluminum cover and front portion covered with clear acrylic canopy to observe safe and dust-free honing.

It is supplied completely with glass honing plate, coarse & fine abrasives and dust cover.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003