

Product Name :
Automatic Slide Projector

Product Code :
NLAB-SCIENTIFICEQ15002



Description :

Automatic Slide Projector

Technical Specification :

Automatic Slide Projector has following specifications:-

- 85 mm thickness.
- 2.5 lenses equipped with Halogen lamp.
- Low pressure cooling system.
- Operated through remote control for forward, reverse and fine focusing.
- 36 slides linear magazine with imported mechanical & optical components to achieve brilliant and sharp projections.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003